


## *Pulse Oximeter Sensors (International)*

*Compatible sensors, adapters, and extension cables that are of equal comfort and quality as those of leading original equipment manufacturers.*



### **Reusable SpO<sub>2</sub> Finger Sensor**

Reusable finger clip sensor is suitable for multiple patient use. It is designed to fit a wide range of digits including adults and pediatrics.



### **Disposable SpO<sub>2</sub> Finger Sensor**

Made from microfoam tape, these adhesive sensors come in four convenient sizes: adult, adult/pediatric, infant, neonate.



### **Soft Sensors**


Soft sensors are designed for multiple patient use for either spot checking or long term monitoring.

Two sizes accommodate a large range of patients from large adults to children as young as age two.



### **Ear Lobe Sensor**

The Earlobe sensor is intended on use with patients who have poor peripheral circulation or when it is not convenient to use a finger sensor. It is especially useful during treadmill tests or when patient motion is anticipated.



### **SpO<sub>2</sub> Y Sensor**

Reusable compatible sensor for use in a variety of SpO<sub>2</sub> monitoring applications.



### **Adapter & Extension Cable**

Allows for greater versatility and convenience with your pulse oximetry sensor.

### **Partial List of OEM Compatibility**

Baxter®, BCI®, Citikon®, Corometrics®, Criticare®, Datascope®, Datex-Ohmeda®, Hewlett Packard®, MEK®, Nellcor®, Nonin®, Novamatrix®, Schiller®, Siemens®, Spacelabs®

**17517 Fabrica Way, Suite H; Cerritos, CA 90703 USA**

**p 714.367.2848 f 714.367.2852**

**[www.mediainc.com](http://www.mediainc.com) [info@mediainc.com](mailto:info@mediainc.com)**