

# Compatible

## Pulse Oximeter Sensors (International)

Compatible sensors, adapters, and extension cables  
 from a wide range of manufacturers.  
 manufacturer.



- **Reusable SpO<sub>2</sub> Finger Sensor**  
 Reusable finger clip sensor is suitable for multiple patient use. It is designed to fit a wide range of digits including adults and pediatrics.




- **Disposable SpO<sub>2</sub> Finger Sensor**  
 Made from microfoam tape, these adhesive sensors come in four convenient sizes: adult, adult/pediatric, infant, neonate.




- **Soft Sensors**  
 Soft sensors are designed for multiple patient use for either spot checking or long term monitoring. Two sizes accommodate a large range of patients from large adults to children as young as age two.



- **Ear Lobe Sensor**  
 The Earlobe sensor is intended on use with patients who have poor peripheral circulation or when it is not convenient to use a finger sensor. It is especially useful during treadmill tests or when patient motion is anticipated.



- **SpO<sub>2</sub> Y Sensor**  
 Reusable compatible sensor for use in a variety of SpO<sub>2</sub> monitoring applications.



- **Adapter & Extension Cable**  
 Allows for greater versatility and convenience with your pulse oximetry sensor.

### Partial List of OEM Compatibility

Baxter®, BCI®, Citikon®, Corometrics®, Criticare®, Datascope®, Datex-Ohmeda®, Hewlett Packard®, MEK®, Nellcor®, Nonin®, Novamatrix®, Schiller®, Siemens®, Spacelabs®  
 Nellcor® is a registered trademark of Covidien LP

17517 Fabrica Way, Suite H; Cerritos, CA 90703 USA

p 714.367.2848 f 714.367.2852

[www.mediaidinc.com](http://www.mediaidinc.com) [info@mediaidinc.com](mailto:info@mediaidinc.com)