

Compatible

Pulse Oximeter Sensors

Compatible sensors, adapters, and extension cables that can be used with the equipment of the original manufacturer.



Reusable SpO₂ Finger Sensor

Reusable finger clip sensor is suitable for multiple patient use. It is designed to fit a wide range of digits including adults and pediatrics.




Disposable SpO₂ Adult Finger Sensor

Made from microfoam tape, these adhesive sensors come in two convenient sizes: adult and adult/pediatric.



Soft Sensors

Soft sensors are designed for multiple patient use for either spot checking or long term monitoring. Two sizes accommodate a large range of patients from large adults to children as young as age two.



SpO₂ Y Sensor

Reusable compatible sensor for use in a variety of SpO₂ monitoring applications.



Adapter & Extension Cable

Allows for greater versatility and convenience with your pulse oximetry sensor.

Partial List of OEM Compatibility

BCI®, Datex-Ohmeda®, Nellcor®, Nonin®
Nellcor® is a registered trademark of Covidien LP

17517 Fabrica Way, Suite H; Cerritos, CA 90703 USA

p 714.367.2848 f 714.367.2852

www.mediaidinc.com info@mediaidinc.com