

MEDIAID

Compatible Veterinary Pulse Oximeter Sensors

Lingual Clip Sensor

The Lingual Clip Sensor is a reusable veterinary sensor suitable for all dogs, cats, and various large animals. The Replaceable silicon pads on this sensor offer added grip to the tongue, which limits sensor movement and reduces interruptions during monitoring. The sensor is packaged with six extra silicon pads. The lingual clip sensor and silicon pads are water resistant and can be easily cleaned.

Rectal Sensor

This veterinary rectal sensor provides pulse oximetry. It is suitable for both large and small animals.

Silicon Pads

The Silicon Pads for the lingual clip are removable pads that can be used multiple times. These pads provide additional adherence to the tongue to limit interruptions while monitoring. The package contains 24 water resistant pads.

OEM Compatibility

Nellcor®, Nonin®, Surgivet®

17517 Fabrica Way, Suite H; Cerritos, CA 90703 USA

p 714.367.2848 f 714.367.2852

www.mediaidinc.com info@mediaidinc.com