


MEDIAID

MediChek

Digital Thermometers



The perfect small digital thermometer for:

Hospitals
Clinics
Physician's Office
Home Health Care
Personal Use

Four Models To Choose From

MediChek Regular – Style Number 6510

Normal response, dual scale with memory and beeper. Ivory battery cap. 5 probe covers, instructions sheets included. Unprinted clear storage case and cap.

MediChek Basal – Style Number 6011

Fast response, dual scale, with flexible probe, memory and beeper. Mauve battery cap. Instruction sheet and 6 graphs included. Opaque storage case and cap.

MediChek Deluxe – Style Number 6012

Fast response, dual scale with flexible probe memory and beeper. Blue battery cap. 5 probe covers and instruction sheet. Opaque storage case and cap.

MediChek Baby – Style Number 6013

Fast response, dual scale with a soft stepped flexible probe, memory and beeper. Red battery cap. Instruction sheet included. Opaque storage case and cap.

MediChek Probe Covers - Style Number 6014

Disposable probe covers minimize the risk of cross contamination.

Equipment Specifications

GENERAL

Display Resolution

1/10 of a degree for both °C & °F

Measurement Range

89.6 °F to 109.4 (32.0 °C to 43.0 °C)

Accuracy

±0.3 °F / ±0.2 °C (89.6 °F to <95.0 °F / 32 °C to <35 °C)
±0.2 °F / ±0.1 °C (95.0 °F to 104.0 °F / 35 °C to 40 °C)
±0.3 °F / ±0.1 °C (>104.0 °F to 109.0 °F / 40 °C to 43 °C)

Measurement Method

Peak Hold System

POWER SOURCE

Battery Type

One 1.5V DC alkaline manganese button size LR41 type or equivalent

Battery Life

Up to 3 years under normal use

WARRANTY

Lifetime Guarantee

17517 Fabrica Way, Suite H; Cerritos, CA 90703 USA

p 714.367.2848 f 714.367.2852

www.mediainc.com info@mediainc.com