

MEDIAID

MediChek

Veterinary Digital Thermometers



● **MediChek Small
Veterinary
Thermometer
- THR010-6013V**

Small flexible insertion tip
- ideal for puppies and kittens!



● **MediChek Large
Veterinary
Thermometer
- THR010-6012V**

Ideal for large animals!

Features

Perfect Utility – For veterinary hospitals, clinics and veterinarian offices

Clinical Accuracy – Fahrenheit and Celsius

Fast Response – In less than 15 seconds

Self-Testing Display – For accurate readings

Visual and Audible Indication – Two ways to let you know when the reading is completed

Specifications

GENERAL

Display Resolution
1/10 of a degree for both °C & °F

Measurement Range
89.6 °F to 109.4 (32.0 °C to 43.0 °C)

Accuracy
±0.3 °F / ±0.2 °C (89.6°F to <95.0°F / 32°C to <35°C)
±0.2 °F / ±0.1 °C (95.0°F to 104.0 °F / 35°C to 40°C)
±0.3 °F / ±0.2 °C (>104.0°F to 109.4°F / >40°C to 43°C)

Measurement Method
Peak Hold System

POWER SOURCE

Battery Type
One 1.5V DC alkaline manganese button size LR41 type or equivalent

Battery Life
Up to 3 years under normal use

WARRANTY
Lifetime Guarantee

17517 Fabrica Way, Suite H; Cerritos, CA 90703 USA

p 714.367.2848 f 714.367.2852

www.mediaidinc.com info@mediaidinc.com