

CROSS REFERENCE CHART FOR VETERINARY COMPATIBLE SENSORS



Compatibility	Monitors	OEM Part Number	Lingual Sensor	Rectal Sensor	Cable Length in Feet	Connector Type		
								
BIONET	BMS VETBM3 VetNext, BM3 VetPlus, OXY9 Vet	B-SPO2-V	CST090-51809	CST080-51809	9 feet (3 meters)	DB9		
CRITICARE	503 Spot, 504 Plus, 504DX, 506DXN, 501, 503, 504, 507, 508, 581, 1100, Poet (602-3A), 8100, Scholar	Not available	CST090-51079	CST080-51079	9 feet (3 meters)	DB9		
DELTALIFE with BCI Technology	OXIPET PLUS - DL410 VET, OXIMETRO OXIPET PLUS DL430 VET, MULTIPARÂMETROS MULT PET - DL900 VET, OXIMETRO OXIPET - DL400 VET	Don't have a specific model, only specifications in Deltalife model	CST090-51309	CST080-51309	9 feet (3 meters)	DB9		
DIXTAL with Novamatrix Technology	DX-2010, DX-2022		CST090-52009	CST080-52009	9 feet (3 meters)	HYPERTRONICS 7 PIN		
EDAN with BCI Technology	VE-H100B	VE-SH7, MS3R-110498	CST090-51309	CST080-51309	9 feet (3 meters)	DB9		
EMAI with BCI Technology	OX-P-10, PX-200, MX-300		CST090-51309	CST080-51309	9 feet (3 meters)	DB9		

CROSS REFERENCE CHART FOR VETERINARY COMPATIBLE SENSORS



Compatibility	Monitors	OEM Part Number	Lingual Sensor	Rectal Sensor	Cable Length in Feet	Connector Type		
								
GOLDWAY with BCI Technology	VET-400A, VET-600A, VET 420F, VET 400C	3078 (BCI Small Lingual Sensor), V-3222 (BCI Large Lingual Sensor)	 CST090-51309	 CST080-51309	9 feet (3 meters)	DB9		
GOLDWAY with Nellcor Technology	VET-400A, VET-600A, VET 420F, VET 400C	VSAT (Lingual Clip Sensor)	 CST090-51209	 CST080-51209	9 feet (3 meters)	DB9		
GUOTENGSCIENCE & TECHNOLOGY DEVELOPMENT Co., Ltd.	GT9003HVET, GT9003MVET		 CST090-52309	 CST080-52309	9 feet (3 meters)	DB9		
HESKA	4400, 4402, 4404, 4700	CAT4060, CAT0451	 CST090-51609	 CST080-51609	9 feet (3 meters)	DB9		
IMFtec with Nellcor Technology	OM-100		 CST090-51209	 CST080-51209	9 feet (3 meters)	DB9		
INSTRAMED with BCI Technology	INMAX VET		 CST090-51309	 CST080-51309	9 feet (3 meters)	DB9		










CROSS REFERENCE CHART FOR VETERINARY COMPATIBLE SENSORS



Compatibility	Monitors	OEM Part Number	Lingual Sensor	Rectal Sensor	Cable Length in Feet	Connector Type		
								
MINDRAY with Mindray Technology	MEC-1200Vet, PM-9000Vet, PM-60Vet, PM-50Vet	0010-10-43237	CST090-51409	CST080-51409	9 feet (3 meters)	REDEL 6 PIN40 DEG DOUBLENOTCH		
MINDRAY with Nellcor Technology	MEC-1000, PM-9000	Not available	CST090-51209	CST080-51209	9 feet (3 meters)	DB9		
MORIYA with BCI Technology	VET 1000		CST090-51309	CST080-51309	9 feet (3 meters)	DB9		
NELLCOR	N-20PAV, N-65VSAT, N-560VSAT, NPV-70, N-85	VSAT (Lingual Clip Sensor)	CST090-51209	CST080-51209	9 feet (3 meters)	DB9		
NEWMED with BCI Technology	OXYLINE		CST090-51309	CST080-51309	9 feet (3 meters)	DB9		
NONIN with Nellcor Technology	9847V, 8500AV, LifeSenseVet, RespSenseVet, PulseSenseVet	2000SL (Lingual Clip Sensor)	CST090-51209	CST080-51209	9 feet (3 meters)	DB9		

CROSS REFERENCE CHART FOR VETERINARY COMPATIBLE SENSORS



Compatibility	Monitors	OEM Part Number	Lingual Sensor	Rectal Sensor	Cable Length in Feet	Connector Type		
								
SURGIVET	V3395, V3402, V3304, V3404, V3204P, V60041, V60046, V90041, V90043, V9200, V9201, V9203, V9204 WWW9230	V1700 (Reflectance Sensor) V1703 (Universal "Y" Lingual Sensor)	CST090-51049	CST080-51049	9 feet (3 meters)	DB9		
Transforma Tecnologia de Ponta Ltda./ECAFIX /FUNBEC with Mediaid Technology	MONITOR CARDIOVET E/ES, OXIMETRO CARDIOVET P, MONITOR DE OXIMETRIA CARDIOVET M	273, 276, 312	CST090-52909	CST080-52909	9 feet (3 meters)	Compushield		
Transforma Tecnologia de Ponta Ltda./ECAFIX /FUNBEC with Mediaid Technology	MONITOR CARDIOVET E/ES, OXIMETRO CARDIOVET P, MONITOR DE OXIMETRIA CARDIOVET M	273, 276, 312	CST090-51909	CST080-51909	9 feet (3 meters)	RJ 12	