


USA & International Cross Reference Chart for Adapter Cables & Extension Cables with Mediaid Part Numbers



Compatibility	Monitors	in 9 Feet Cable Length	Connector Type			
			Monitor side		Sensor side	
BCI	3100, 3101, 6100, 9100, Advisor, Autocorr, Minitorr, Minitorr Plus	CST110-4101	REDEL 7 PIN SINGLE NOTCH		DB9	
BCI	3300, 3301, 3302, 3303, 3304, 3401	CST110-4102	DB9		DB9	
CRITICARE	503 Spot, 504Plus, 504DX, 506DXN	CST110-4105	DB9		DB9	
CRITICARE	501, 503, 504, 507, 508, 581, 1100, Poet (602-3A), 8100, Scholar	CST110-4106	LEMO 5 PIN		DB9	
DATASCOPE	Accusat, Accutorr 3/4 SAT, Accutorr Plus (Datascope SpO2), Passport (Datascope SpO2)	CST110-4110	DIN 8		DB9	
DATEX	AS/3, CS/3, CH/S, CH/RS, CG/2S, Cardiocap II, Capnomac, Ultima, Light Monitor, M-Estpr, Oscar II, OSP-200, Satelite Trans Plus, ULTSV I, ULTSVO, ULTSIO, ULTS00	CST110-4111	DPK1		DB9	
HEWLETT PACKARD	78352C, 78354C, 78834C, M1020A, M1025B, M1025A, M1350B, M1722A/B, M1732A/B, M2475B, Virida A1, A3	CST110-4113	DPK1		DB9	
HEWLETT PACKARD	78352C, 78354C, 78834C, M1020A, M1025B, M1025A, M1350B, M1722A/B, M1732A/B, M2475B, Virida A1, A3	CST110-4133	Amphinol		DB9	

USA & International Cross Reference Chart for Adapter Cables & Extension Cables with Mediaid Part Numbers



Compatibility	Monitors	in 9 Feet Cable Length	in 9 Feet Cable Length			
			Monitor side		Sensor side	
MARQUETTE EAGLE	Eagle, Dash, SOLAR, TRAM Modular 250SL, 450SL, 650SL, 850A, 850SL	CST110-4116	MARQUETTE Flat 11 PIN		DB9	
MARQUETTE HELLIGE	Hellige Eagle 1000, 3000, SMU, SMK, SMV	CST110-4117	HELLIGE ROUND 7 PIN		DB9	
MARQUETTE HELLIGE	7000 Series, TRAM A.R., TRAM 200SL, 400SL, 600SL, 800A, 800SL	CST110-4127	HYPERTRONICS 7 PIN		DB9	
MASIMO	All Models	CST110-4118	3M		DB9	
NELLCOR®	Adapter Nellcor C-Cal to Oximax	CST110-4163	DB9		DB9	
NELLCOR®	N100, 185, 195	CST110-4104	DB9		DB9	
NELLCOR®	NPB-290, 295, N390, 395, 3000, 6000	CST110-4120	3M		DB9	
NELLCOR®	N180, 200, NPB-190, N6000, NPB-4000,3900, 3940, 3910, 75	CST110-4146	DB9		DB9	
NIHON KOHDEN	AL-800PA Module, BSM-1101, 2101/2, 3101, 7100, 7200, OLV-1100A, Lifescope	CST110-4121	DPK1 10 PIN		DB9	

USA & International Cross Reference Chart for Adapter Cables & Extension Cables with Mediaid Part Numbers



Compatibility	Monitors	in 9 Feet Cable Length	in 9 Feet Cable Length			
			Monitor side		Sensor side	
NONIN	8500	CST110-4122	DB9		DB9	
NOVAMETRIX	505, 510, 511, 515B/C, 520, 7000, 7100, Oxypleth	CST110-4124	HYPERTRONICS 6 PIN		DB9	
NOVAMETRIX	505, 510, 511, 515B/C, 520, 7000, 7100, Oxypleth	CST110-4125	HYPERTRONICS 6 PIN		HYPERTRONICS	
OHMEDA	3700, 3710, 3770, 3775, Biox 3740, 3760, 3800, 3900, 4700, Tuffsat, OxiCap, 5250 RGM, Rascal II	CST110-4147	HYPERTRONICS 7 PIN		DB9	
SIEMENS	Sirecust 700, 630, Sirem SpO2 Module	CST110-4130	SIEMENS (BIG ROUND)		DB9	
SIEMENS		CST110-4138	SIEMENS (SMALL ROUND)		DB9	
SPACELABS	90419-01, 90651A-07	CST110-4124	HYPERTRONICS 6 PIN		DB9	
SPACELABS	90351-0/6, 90465, 90466, 90467, PC Express, 90489, 90651A-08, Im77	CST110-4132	HYPERTRONICS 7 PIN		DB9	
MINDRAY	MEC-1000, PM-9000	CST110-4152	REDEL 6 PIN 40 DEG DOUBLE NOTCH		DB9	